

Supplemental Table 2: Primers for ChIP assay.

S.No	Name	Primers
1	<i>Myc</i>	Sense: 5'-GCC CTT TCC CCA GCC TTA GC-3' Antisense: 5'-AAC CGC ATC CTT GTC CTG TGA GTA-3'
2	<i>miR-23a/27a</i>	Sense: 5'-TAG AGGACC ACG TGT TCA TTT TGC-3' Antisense: 5'-TGT AAA ATG GAGCCA AGA GCC TC-3'
3	<i>miR-17/92</i>	Sense: 5'-ACT TTG CAG TCT CGG GTG TTC-3' Antisense: 5'-CGG GAT AAA GAG TTG TTT CTC CAA-3'